

Account Name: \_\_\_\_\_ Patient Name: \_\_\_\_\_  
Doctor: \_\_\_\_\_ Weight: \_\_\_\_\_ Shoe Size: \_\_\_\_\_ Age: \_\_\_\_\_  
Address: \_\_\_\_\_ Gender: \_\_\_\_\_  
\_\_\_\_\_ Diagnosis: \_\_\_\_\_  
City: \_\_\_\_\_  
State: \_\_\_\_\_ Zip: \_\_\_\_\_  
Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ Casted By: \_\_\_\_\_ Date: \_\_\_\_\_



- Type:**  Right  Left  Bilateral
- Color:**  
 Black  White  Chocolate Brown  Light Beige  
 Tan  Pink  Dark Beige (Taupe)  Powder Blue
- Closure:**  Lace Only  Velcro Only  Combination (Laces with Velcro Strap)
- Brace Height:**  6" Rolled Collar  7"  9"  As Marked On Cast
- Reinforcement:**  Flexible (2mm)  Semi-Rigid (3mm)  Firm (4mm)
- Heel:**  
 Standard (Solid Plastic Heel)  No Reinforcement (Leather Covered Heel, No Plastic)
- Ankle Corrections:**  
 Correct Anterior/Posterior to 90°  As Casted  
 Correct Medial/Lateral to 90°  As Casted
- Forefoot Corrections:**  Correct to 90°  As Casted
- Poly Footplate Length:**  
 Proximal to Met Heads  Sulcus  Full Foot
- Fabric Footplate Length:**  
 Proximal to Met Heads  Sulcus  Full Foot

Special Cast/Brace Modification Instructions: \_\_\_\_\_

Casting Material (STS Mid-Leg Casting Sock. Additional charge of \$13.00 each. Please specify size and quantity):

None  \_\_\_\_\_ Small  \_\_\_\_\_ Medium  \_\_\_\_\_ Large  \_\_\_\_\_ X-Large

**Casting Reminders:**

- Cast should be taken with the ankle at 90° and the foot in subtalar neutral.
- Patient name and Doctor/practice name must be written on the cast.
- Marks should be made on the cast on any areas that need special attention and noted on this order form on how you would like them to be accommodated.

**Additional Charges:**

Corrections or changes to the brace after fabrication may incur additional charges. To avoid these charges please review your cast and order form for accuracy before shipping.

Toll Free: 866-406-7926  
Email: [rybomedical@gmail.com](mailto:rybomedical@gmail.com)  
Production Facility:  
175 Park Street  
Livermore Falls, ME 04254

Fax: 949-446-8224  
[www.rybomedical.com](http://www.rybomedical.com)  
Corporate Office:  
21302 Calle Horizonte  
Lake Forest, CA 92630